
Puralle





Class II, Type A2 Biological Safety Cabinets

- Creating Safe Solutions for Life Science Laboratories
- Engineered for Simplicity and Efficiency
- Certified for Safety and Performance



230 watt1 Purair BIO model AS-AHA-133-CA-A



220-490 watt1

The single EC blower motor assures lower cost of ownership in one of the world's most energy efficient biological safety cabinets.





NSF Certification applies to the biological safety cabinet models AHA-133-AA-B (Air Science Model (AS-AHA-133-AA-B), AHA-133-AB-B) are Science Model AS-AHA-133-AB-B) and AHA-133-CB-B (Air Science Model AS-AHA-133-CB-B) and AHA-133-CB-B) and AHA-133-CB-B) and AHA-133-CB-B) manufactured by Chung Fu Technical Development Co. Ltd., Taipei, Taiwan and marketed by Air Science, LLC, USA.





INTRODUCTION

The Purair® BIO biological safety cabinet (BSC) is designed for safety and performance in accordance with US and International Standards such as NSF/ANSI 49 which certify that Class II, Type A2 laminar flow cabinets are suitable for working with biosafety agents at levels 1, 2, 3 and 4. Air Science 4-foot models are certified by NSF.

HEPA filtration of downflow and exhaust paths provides a primary containment work area for life science research, cell culture processing and other applications where protection of the user, the work product, and the environment and mitigation of cross-contamination on the work surface are needed.



APPLICATIONS

Purair BIO is designed to protect individuals, the environment and products from a variety of biological particulates. Specific applications include, but are not limited to:

Life Science Research \ Sterile Product Preparation \ Biological Protocols





- The Purair BIO does not use a costly and overcomplicated microprocessor-based controller. Our basic electronic control system provides simple and reliable oversight of all cabinet systems.
- Our design is based on a single, energy-efficient EC Brushless DC Motor and air circulation system which manages airflow at all critical points including inflow, downflow and exhaust. There is no need for dual motor synchronization, balancing and expense.
- Front access to both HEPA supply and exhaust filters encourages quick, safe filter removal and replacement by an authorized technician.

ENHANCED PROTECTION, EASE OF USE

The Purair BIO maintains negative pressure inside the cabinet during operation to prevent contaminants from escaping the work area. HEPA filtration scours 70% of the incoming room air to protect the products, while the remaining 30% of the exhausted air is filtered by a second HEPA filter. Purair BIO cabinets provide ample workspace with environmentally sound operation, low energy consumption and user-friendly operation.

Single EC Blower Motor Design. The Purair BIO is designed with a single EC blower motor for ease of use, reliability and to promote low cost of ownership. Dual motor designs can introduce a number of problems, including unbalanced airflow, higher maintenance costs, longer certification processes and uneven filter loading. Airflow adjustments are simpler, operational costs lower and maintenance easier for BSCs equipped with single EC blower motors.

Flexibility. The Purair BIO includes multiple service connections for maximum flexibility. It includes duplex electrical outlets with splash-proof, UL listed covers. CSA certified service valves for gas, air and vacuum are also available. Maximum working pressure of 75PSI.

Operator Safety. The anti-ultraviolet, 6 mm tempered glass ensures maximum protection for the operator.

Ergonomic Design. The user interface combines ergonomics, safety and aesthetics with a 10° angled window design that reduces operator head and elbow discomfort, as well as eye strain and fatigue.



model AS-AHA-133-CA-B.

The slim profile design passes easily through a standard 38" wide doorway to accommodate a variety of new or replacement installations.



This product exceeds OSHA, ANSI and other International Certification Standards. Specifications are subject to change without notice or obligation on the part of Air Science. For questions contact Air Science.

¹⁾ Energy consumption disclosure is based on internal testing with the blower running and lights on.

Design Features (p.3) Performance & Selection (p.4) Containment & Filtration (p.5)



DESIGN FEATURES

- **A.** Energy Efficient: The quiet, internal EC Brushless
- **B.** Standard Control System: Includes simple, reliable
- C. HEPA Filtration Lock: The patented Quick Access
- **D.** Paper Catch: Protective screen located at the
- **E.** Double-Wall Plenum Design: Double-wall design creates a unique plenum which surrounds contami-
- F. Air Velocity: The air velocity and associated
- **G.** Surround Air Intake Grille: All contaminated air is
- H. Outlet: Outlets are installed in the work area to
- I. Durable Interior: The Purair BIO utilizes a heavy-

- **J.** Ergonomic Fit: The angled front, narrow-front
- K. Safe: Includes HEPA filters (Class H14) tested to a

ADDITIONAL FEATURES

Robust Cabinet Construction: Key components, including

NSF 49 & EN12469 Listed: NSF 49 (4-foot models only)

Safety Interlock Feature: Activates the fan and fluorescent

Standards Compliant: The Purair BIO undergoes 10

Full Product Support: The best value Class II Biological

220 watt¹ Purair Bio model AS-AHA-103-CA-B.

This product exceeds OSHA, ANSI and other International Certification Standards. Specifications are subject to change without notice or obligation on the part of Air Science. For questions contact Air Science.

¹⁾ Energy consumption disclosure is based on internal testing with the blower running and lights on.

The Purair BIO includes industry-leading innovations and technology. It is easy to install, energy efficient, cost effective, and safe.

PERFORMANCE

The supply/exhaust HEPA filters provide 99.99% efficiency at 0.3 microns (Class H14).

DESIGN


The Purair BIO is self-contained and does not require venting to the outside. Because filtered air is returned to the room, there is no increased load requirement for facility HVAC make-up air. This eliminates the cost of additional facility ductwork construction as well as HVAC maintenance and overhead.

RFI IABILITY

An innovative surrounded air intake grille keeps all contaminated air contained within the cabinet and isolated thanks to a double-wall design and negative internal pressure.



PERFORMANCE & SELECTION



SELECTION

Purair BIO biological safety cabinets are available in 4 standard sizes.

The slim profile hood design is the thinnest on the market and allows multiple laboratory configurations. Units are portable and may be moved from one location to the next with minimal downtime and without filter changes. Set-up, operation and filter maintenance are straightforward.


CONTROL

Purair BIO biological safety cabinet uses time-proven electronic controls instead of costly microprocessor controllers.

Because all Class II, Type A2 cabinets must meet NSF standards for airflow, face velocity and other performance attributes, eliminating superfluous control and indicating devices simplifies operation and user interface while essential functions are maintained. As a result, the Purair BIO Series offers a cost-effective, efficient and compliant solution to biological safety cabinet applications in life science and associated uses.

All Purair BIO EN models are equipped with membranestyle switches.

- A. An analog Minihelic[™] gauge verifies airflow.
- B. On/Off key-switch controls main power.
- C. Soft key switches control interior lighting, UV lamp, interior outlet activation and fan.
- D. Sash alarm indicator activates an audible and visual warning of unsafe sash position (10 inches or 25 cm from the bottom of cabinet).
- E. Mute button silences the alarm for a specified time.



This product exceeds OSHA, ANSI and other International Certification Standards. Specifications are subject to change without notice or obligation on the part of Air Science. For questions contact Air Science.

¹⁾ Energy consumption disclosure is based on internal testing with the blower running and lights on

CONTAINMENT AND PROTECTION

The Purair BIO maintains an airflow ratio of 70% recirculation to 30% exhaust to ensure operator protection. The inflow and downflow balance is precisely established, with no room air entering the work zone to prevent product contamination. Airflow patterns are precisely tuned and tested to create an optimum air curtain on the front aperture, maintaining personnel and product protection even during the unlikely event of severe inflow or downflow imbalances.

The integrated HEPA filtration system provides clean air to the work surface in a gentle vertical laminar flow pattern, allowing the exhaust HEPA filter to trap biohazardous particles prior to the air being exhausted into the room.



What We Avoid

- Costly microprocessor controllers
- Dual blower motors with twice the energy consumption



HEPA FILTRATION

The Purair BIO uses HEPA filters to provide a range of high performance protection.

These self-contained filters are designed to physically capture particles larger than 0.3 microns with >99.995% typical efficiency.

The filters feature an integral groove filled with gel at the air inlet side, ensuring a perfect seal to the housing system. The aluminum frame guards against swelling typical of wooden framed filters.

For unobstructed airflow and superior filtration, filters do not contain aluminum separators.

A patented HEPA filtration lock maintains filter efficiency, minimizes the chance of leakage and prolongs filter life. The filters can also be changed from the front side of the cabinet quickly and easily.

ENERGY EFFICIENCY

The Purair BIO maintains one of the world's highest performance ratings for a brushless DC motor. Additional benefits include:

Better Apportioned Power. Over 80% of the EC motor output power is converted to kinetic energy to ensure sustainable energy savings over the life of the motor.

Extended Filter Life. Balanced airflow and even distribution of downflow and exhaust paths promotes uniform filter loading to prolong filter life.


Constant Feedback Motor Speed. The EC motor automatically adjusts speed to maintain compliant airflow at all critical points while compensating for filter loading and facility voltage fluctuations.

AIRFIOW

The Purair BIO is configured to comply with either NSF 49 or EN12469 criteria for airflow within critical points of the cabinet. In either model, the combination of HEPA supply and exhaust filters yields a fully integrated performance envelope for product, personnel and environmental protection from particulates.

Face Velocity (Cabinet Intake)	NSF 49	EN12469	
Note: 105 fpm is the midpoint	~100 fpm	~80 fpm	
for the approved range.	~0.5 m/s	~0.4 m/s	

An optional exhaust collar may be added which allows the Purair BIO to be connected to a facility exhaust ventilation system.



This product exceeds OSHA, ANSI and other International <u>Certification</u> Standards.

Specifications are subject to change without notice or obligation on the part of Air Science. For questions contact Air Science.

¹⁾ Energy consumption disclosure is based on internal testing with the blower running and lights on.

6

Specifications (p.6)
Options & Accessories (p.8)

AS-AHA-103

3 ft.

AS-AHA-133




AS-AHA-163



AS-AHA-193





MODEL	VOLTAGE	CE	RTIFICAT	IONS	PROTECTION		DIMENSIONS			WEIGHT (LBS/KG)	
		NSF 49	EN12469	IEC61010-1	Particulates* (Personnel, Work Surface, Environment)	Vapors (Personnel, Work Surface)	Internal (W × D × H)	External (W × D × H)	Shipping (W × D × H)	Net	Ship

Purair BIO Models

Furair BIO Models											
3 ft. AS-AHA-103-CA-B 3 ft. AS-AHA-103-CB-B	115V, AC, 60Hz 230V, AC, 50Hz	No	Yes	Yes	Yes	When ducted to outside	36.2" × 24.6" × 28" 920 × 626 × 709 mm	40.7" × 31.1" × 60.5" 1034 × 789 × 1537 mm	43.3" × 33.5" × 77.8" 1100 × 850 × 1950 mm	463 / 210	529 / 240
4 ft. AS-AHA-133-CA-A 4 ft. AS-AHA-133-CB-A	115V, AC, 60Hz 230V, AC, 50Hz	No	Yes**	Yes	Yes	When ducted to outside	48" × 24.6" × 28" 1220 × 626 × 709 mm	52.5" × 31.1" × 60.5" 1334 × 789 × 1537 mm	55.1" × 33.5" × 77.8" 1400 × 850 × 1950 mm	529 / 240	617 / 280
4 ft. AS-AHA-133-CA-B 4 ft. AS-AHA-133-CB-B	115V, AC, 60Hz 230V, AC, 50Hz	Yes	No	Yes	Yes	When ducted to outside	48" × 24.6" × 28" 1220 × 626 × 709 mm	52.5" × 31.1" × 60.5" 1334 × 789 × 1537 mm	55.1" × 33.5" × 77.8" 1400 × 850 × 1950 mm	529 / 240	617 / 280
5 ft. AS-AHA-163-CA-B 5 ft. AS-AHA-163-CB-B	115V, AC, 60Hz 230V, AC, 50Hz	No	Yes	Yes	Yes	When ducted to outside	59.8" × 24.6" × 28" 1520 × 626 × 709 mm	64.3"x 31.1" × 60.5" 1634 × 789 × 1537 mm	67" × 33.5" × 77.8" 1700 × 850 × 1950 mm	617 / 280	705 / 320
6 ft. AS-AHA-193-CA-B 6 ft. AS-AHA-193-CB-B	115V, AC, 60Hz 230V, AC, 50Hz	No	Yes	Yes	Yes	When ducted to outside	71.6" × 24.6" × 28" 1820 × 626 × 709 mm	76.1" × 31.1" × 60.5" 1934 × 789 × 1537 mm	78.7" × 33.5" × 77.8" 2000 × 850 × 1950 mm	705 / 320	816 / 370

^{* 99.99% @ 0.3} microns

^{** 4} ft. EN versions have a 20 cm window position, all other models have a 25 cm window position.

(2) T5, 35 watt

(1) FL 40 watt, 254 nm

< 60

(4) T5, 21 watt

(2) FL 20 watt, 254 nm

< 62

AS-AHA-133-CA-A AS-AHA-103-CA-B AS-AHA-133-CA-B AS-AHA-163-CA-B AS-AHA-193-CA-B **Filtration** AS-AHA-103-CB-B AS-AHA-133-CB-B AS-AHA-133-CB-A AS-AHA-163-CB-B AS-AHA-193-CB-B Airflow <... 70% recirculated / 30% exhausted ...> Average Inflow Velocity <... 100 fpm (0.51 m/s) ...> Average Downflow Velocity 64 fpm (0.325 m/s) 60 fpm (0.30 m/s) 60 fpm (0.30 m/s) 60 fpm (0.30 m/s) 60 fpm (0.30 m/s) Exhaust Air Volume with Exhaust Canopy 353 cfm (10 m³) 459 cfm (13 m³) 388 cfm (11 m³) 565 cfm (16 m³) 671 cfm (19 m³) Supply HEPA Filter <... (1) 99.99% @0.3 microns (Class H14) ...> Exhaust HEPA Filter <... (1) 99.99% @0.3 microns (Class H14) ...> AS-AHA-103-CA-B AS-AHA-133-CA-B AS-AHA-133-CA-A AS-AHA-163-CA-B AS-AHA-193-CA-B Construction AS-AHA-103-CB-B AS-AHA-133-CB-B AS-AHA-133-CB-A AS-AHA-163-CB-B AS-AHA-193-CB-B Room Height 89.4" / 227 cm Construction <... SUS#304; 6 mm tempered window glass ...> Window Operation 9.8" / 25 cm 9.8" / 25 cm 7.9" / 20 cm 9.8" / 25 cm 9.8" / 25 cm Blower <... EC brushless DC Motor ...> **Electrical Switches** <... Membrane push switches; On/Off ...> **Electrical Outlets** (2) Electrical outlets (2) Duplex electrical outlets with (2) Electrical outlets (2) Electrical outlets (2) Electrical outlets (CE certified) splash-proof cover (UL listed) (CE certified) (CE certified) (CE certified) Electrical <... 115V, 60Hz, 15A or 230V, 50Hz, 10A ...> AS-AHA-103-CA-B AS-AHA-133-CA-B AS-AHA-133-CA-A AS-AHA-163-CA-B AS-AHA-193-CA-B Efficiency AS-AHA-103-CB-B AS-AHA-133-CB-B AS-AHA-133-CB-A AS-AHA-163-CB-B AS-AHA-193-CB-B Power Consumption 220 watt 300 watt 230 watt 400 watt 490 watt

PRODUCT SPECIFICATIONS

Fluorescent Intensity lux

Fluorescent Lamp

UV Lamp

Noise, dBA

<... > 1,000 ...>

(2) T5, 28 watt

(1) FL 30 watt, 254 nm

< 57

(2) T5, 28 watt

(1) FL 30 watt, 254 nm

< 63

(2) T5, 21 watt

(1) FL 20 watt

< 58

Specifications are subject to change without notice.

OPTIONS & ACCESSORIES

Purair BIO Models		AS-AHA-103-CA-B AS-AHA-103-CB-B	AS-AHA-133-CA-B AS-AHA-133-CB-B AS-AHA-133-CA-A AS-AHA-133-CB-A	AS-AHA-163-CA-B AS-AHA-163-CB-B	AS-AHA-193-CA-B AS-AHA-193-CB-B
Exhaust Collar	For connecting the biosafety cabinet to building exhaust or remote blower for outside ducting.	HBE-103-AA-A	HBE-133-AA-A	HBE-163-AA-A	HBE-193-AA-A
Service Fixture Kit	Additional connections available for Air/Gas/Vacuum. CSA certified.	BA-L4100-158-B	BA-L4100-158-B	BA-L4100-158-B	BA-L4100-158-B
IV Bar Kit	Includes 6 hooks. Specify when ordering. Field installed.	BC-1934-A	BC-1671-A	BC-1673-A	BC-1674-A
Moveable Elbow / Arm Rest	Made of stainless steel #304 and provides ergonomic forearm support to prevent grill blockages and improve posture.	BC-1933-A	BC-1933-A	BC-1933-A	BC-1933-A

ELECTRICAL SAFETY AND CERTIFICATION

- All components meet or exceed applicable safety requirements.
- Each cabinet is individually tested for electrical safety at the factory and documentation specific to each tested cabinet is maintained on file.
- IEC Certified for USA and Canada
- Models available with NSF 49 and EN12469 certifications.

Contact Air Science or your sales rep for preparation information.

WARRANTY

The Purair BIO and BIO-EN are supported by a standard industry warranty which starts on the date of shipment from our factory. Under the protection period this warranty covers defects in materials and workmanship. Our liability is, at our option, to repair or replace any defective parts of this equipment if you document that these parts were defective at the time we sold the product to you.

This Warranty is valid on specified Air Science products with regard to mechanical parts AS LONG AS genuine Air Science parts and filters were used in compliance with Air Science specifications.

This Warranty constitutes the customer's sole and exclusive remedy and Air Science's sole and exclusive responsibility with respect to any alleged breach of this limited warranty.

For details visit the Service section of our website at www.airscience.com.

STANDARDS AND COMPLIANCE					
Quality Management Systems	ISO 9001:2015				
Cabinet Performance	IEST-RP-CC002.2 AS 2252.6				
Air Quality	ISO 14644-1, Class 4				
Sterile Compounding	USP 797				
Filtration	IEST-RP-CC034.1 IEST-RP-CC001.3 EN 1822				
Electrical Safety	IEC61010-1				



